

### Product Overview

- Adjustable height from 57 cm to 87 cm
- 2 liftable columns
- Trespa™ chemical resistant worksurface W45 cm x D75 cm
- Weight capacity: 250 kg
- Worksurface can easily be lifted up or down by commuting an electrical switch
- Worldwide electrical compatibility
- Recommended use with liquid chromatograph
- HPLC Ionbench is on caster wheels and can easily be moved
- Convenient access to the top of your HPLC or UHPLC
- Safe access to solvents

Liftable  
motorized  
worksurface



### *Control :*



#### Standard switch


- *Up and Down*



#### Optional position display with 3 memories (ref : BCHPOS)

- *Up and Down*
- *Memory function ( 3 in total )*
- *Display for height indication*

### *Photo example :*



### *Detailed Technical specifications :*

- Compatible with 230V/50HZ or 110V/60HZ
- Electronic components are all compliant with RoHS directives (Restriction of Hazardous Substances Directive).

External Dimensions	BCHLC4575
Depth	75 cm / 29,5 inch
Width	45 cm / 17,7 inch
Height	From 57 to 87 cm / From 22,4 to 34,3 inch
Weight :	30 kg / 66 lb

### *Summary*

Description	Part-number
Ionbench-LC BCHLC4575 / 220 Volt	BCHLC4575V220
Ionbench-LC BCHLC4575 / 110 Volt	BCHLC4575V110
Optional position display with 3 memories	BCHPOS