



OIL DRAIN PUMP FOR VACUUM PUMPS (ODP)

Product Overview

- 1 small compact oil pump with on/off switch and power supply 220V/110V (worldwide compatible plugs provided).
- 2 plastics hoses for oil pumping into a waste container (empty bottle, ...)
- Simply connect the inlet pipe to the MS Noise Oil Drain Kit
- Place the outlet pipe into the waste container
- Turn on the oil drain pump
- When your pump is empty, just turn off the oil pump
- Oil can be changed without any tool



Description	Part-number
Oil drain Pump	ODP
Oil drain Kit	ODK-XXX

MS NOISE EUROPE & ASIA - MS NOISE - Building 4 ZAC des Vauguilletes - 1, bd des Noyers Pompons - F-89100 SENS – France (contact@msnoise.com)
MS NOISE NORTH AMERICA - MS NOISE - Suite 410 - 185, Alewife Brook Parkway - Cambridge, MA 02138 – USA (contact_US@msnoise.com)