












All ESI Racks are Polypropylene or UHMW polyethylene, are highly chemically resistant, and steam autoclavable at 121° C (250° F).

**Large Racks** Dimensions: 247mm L x 101mm W x 74mm H (9.72" x 3.98" x 2.91").  
Ships flat - Easy to assemble.

Rack P/N	Description	Compatible A/S	Vials	Cover
LR-90-13 	90 positions 6 x 15 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC-8, SC14	8mL 13mm dia	N/A
LR-90-13R 	90 positions round holes 6 x 15 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC-8, SC14	8mL 13mm dia	N/A
LR-60-16 	60 positions 5 x 12 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC-8, SC14	12mL, 15mL 16mm dia	N/A
LR-40-20 	40 positions 4 x 10 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC-8, SC14	20mm dia	N/A
LR-21-30 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC-8, SC14	30mL, 50mL, 30mm dia	N/A

## Super Racks Dimensions: 130mm L x 41mm W (5.12" x 3.98").

Rack P/N	Description	Compatible A/S	Vials	Cover
SR2-120-13 	120 positions 6 x 20 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	8mL 13mm dia	N/A
SR4-60-14 	60 positions 5 x 12 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	1mL, 2mL, 4mL 14mm dia	SR4-00-01 
SR2-80-16 	80 positions 5 x 16 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	13mL 16mm dia	N/A
SR2-21-22 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	12mL 22mm dia	N/A
SR2-27-22 	27 positions 3 x 9 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	12mL, 22mm dia	N/A


## Super Racks (continued)

Rack P/N	Description	Compatible A/S	Vials	Cover
SR2-21-27 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	20mL, 25mL 27mm dia	N/A
SR2-27-28 	27 positions 3 x 9 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	20mL, 50mL 28mm dia	N/A
SR2-21-30 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	30mL, 50mL 30mm dia	N/A
SR2-27-31 	27 positions 3 x 9 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	31mm dia	N/A
SR2-21-35 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	35mm dia	N/A







## Super Racks (continued)

Rack P/N	Description	Compatible A/S	Vials	Cover
SR2-12-40 	12 positions 2 x 6 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	40mm dia	N/A
SR2-12-50 	12 positions 2 x 6 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	125mL 50mm dia	N/A
SR14-40-50 	40 positions 4 x 10 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	125mL 50mm dia	N/A
SR2-10-61 	10 positions 2 x 5 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	250mL 61mm dia	N/A
SR14-40-61 	40 positions 4 x 10 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	250mL 61mm dia	N/A







## Super Racks (continued)

Rack P/N	Description	Compatible A/S	Vials	Cover
SR14-21-72 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4	72mm dia	N/A





## Micro Racks Dimensions: 158.7mm L x 67.3mm W (6.25" x 2.65")

Rack P/N	Description	Compatible A/S	Vials	Cover
MR-60-08 	60 positions 5 x 12 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	0.5mL 8mm dia	MR-00-01 
MR-90-08 	90 positions 6 x 15 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	0.5mL 8mm dia	MR-00-01 
MR-21-14 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	1mL, 2mL, 4mL 14mm dia	MR-00-01 

### Micro Racks (continued)











Rack P/N	Description	Compatible A/S	Vials	Cover
MR-40-14 	40 positions 4 x 10 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	1mL, 2mL, 4mL 14mm dia	MR-00-01 
MR-21-15 	21 positions 3 x 7 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	1mL, 2mL, 4mL 15mm dia	MR-00-01 
MR-40-15 	40 positions 4 x 10 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	1mL, 2mL, 4mL 15mm dia	MR-00-01 

### Microtiter Plates Dimensions: 127.4mm L x 85.1mm W (5.01" x 3.35").

Rack P/N	Description	Compatible A/S	Cell Volume	Cover
MT-96-500-05-C 	96 cells 8 x 12 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC- $\mu$	500 $\mu$ L	MT-00-02 
MT-96-2mL-02 	96 cells 8 x 12 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC- $\mu$	2mL	MT-00-02 



## Microtiter Plates (continued)

Rack P/N	Description	Compatible A/S	Cell Volume	Cover
MT-96-4mL-02 	96 positions 8 x 12 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4 SC- $\mu$	4mL	MT-00-02 
MT-48-5mL-02 	48 positions 6 x 8 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	5mL	MT-00-02 
MT-48-7.5mL-02 	48 positions 6 x 8 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	7.5mL	MT-00-02 
MT-48-10mL-02 	48 positions 6 x 8 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	10mL	MT-00-02 
MT-24-10mL-02 	24 positions 4 x 6 matrix  Applications:	SC-2, SC-4, SC-E2, SC-E4, SC- $\mu$	10mL	MT-00-02 

## Standards Racks

Rack P/N	Description	Compatible A/S	Vials	Cover
<p>ST10-28</p> 	<p>10 positions 6+4 matrix</p> <p>Applications:</p>	<p>SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14</p>	<p>25mL, 50mL 28mm dia</p>	<p>N/A</p>
<p>ST10-28-CP</p> 	<p>Compact Spacing 10 positions 6+4 matrix</p> <p>Applications:</p>	<p>SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14</p>	<p>25mL, 50mL 28mm dia</p>	<p>N/A</p>
<p>ST12-14</p> 	<p>12 positions 7+5 matrix</p> <p>Applications:</p>	<p>SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14,</p>	<p>1mL, 2mL, 4mL 14mm dia</p>	<p>N/A</p>
<p>ST12-16</p> 	<p>12 positions 7+5 matrix</p> <p>Applications:</p>	<p>SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14</p>	<p>13mL 16mm dia</p>	<p>N/A</p>
<p>ST10-1</p> 	<p>Corzan 10 positions 6+4 matrix</p> <p>Applications:</p>	<p>SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14</p>	<p>25mL, 50mL 28mm dia</p>	<p>N/A</p>



## Standards Racks (continued)

Rack P/N	Description	Compatible A/S	Vials	Cover
<p>ST10-28-EL</p> 	<p>10 positions 5+5 matrix</p> <p>Applications:</p>	SC-EL	25mL, 50mL 28mm dia	N/A
<p>ST10-28-MR</p> 	<p>10 positions 6+4 matrix</p> <p>Applications:</p>	SC- $\mu$	25mL, 50mL 28mm dia	N/A
<p>ST5-28</p> 	<p>5 positions</p> <p>Applications:</p>	SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14	25mL, 50mL 28mm dia	N/A
<p>ST10-31</p> 	<p>10 positions 6+4 matrix</p> <p>Applications:</p>	SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14	30mL, 50mL 31mm dia	N/A
<p>ST9-34</p> 	<p>9 positions 5+4 matrix</p> <p>Applications:</p>	SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14	34mm dia	N/A

## Standards Racks (continued)

Rack P/N	Description	Compatible A/S	Vials	Cover
ST5-61	5 positions Applications:	SC-2, SC-4, SC-E2, SC-E4, SC-8, SC14	250mL 61mm dia	N/A



Rack P/N	Description	Compatible A/S	Vials	Cover
ST7-60	7 positions Applications:	SC-4	250mL 61mm dia	N/A



2010.09.30