

Product Overview

- Adjustable height from 57 cm to 87 cm
- 1 liftable Column
- Trespa™ chemical resistant worksurface W45 x D55 cm
- Weight capacity: 150 kg
- Worksurface can easily be lifted up or down by commutating an electrical switch
- Worldwide electrical compatibility
- Recommended use with liquid chromatograph
- HPLC Ionbench is on caster wheels and can easily be moved
- Convenient access to the top of your HPLC or UHPLC
- Safe access to solvents

Liftable

motorized

worksurface



Control :



Standard switch


- Up and Down



Optional position display with 3 memories (ref : BCHPOS)

- Up and Down
- Memory function (3 in total)
- Display for height indication

Photo example :



Detailed Technical specifications :

- Compatible with 230V/50HZ or 110V/60HZ
- Electronic components are all compliant with RoHS directives (Restriction of Hazardous Substances Directive).

External Dimensions	cm	inch
Depth	65	25.59
Width	55	21.65
Height	From 57 to 87 cm	From 22,4 to 34,3
Weight :	22 kg	48.50 lb

Summary

Description	Part-number
Ionbench-LC BCHLC4555 / 220 Volt	BCHLC4555V220
Ionbench-LC BCHLC4555 / 110 Volt	BCHLC4555V110
Optional control position display with 3 memories	BCHPOS