P/N ABAPC00X



ECB2004PB

GRADIENT BOX WITH DEGASSER AND PC

Box is suited for liquid handling in ECOM gradient systems. It accommodates a container for solvent bottles, built-in highefficiency on-line degasser and four-way gradient valve. To have HPLC system as small as possible this box also includes industrial PC board. This solution allows installation of OS (Windows or Linux) and chromatography SW ECOMAC or Clarity.

Degasser - unique design assures reliable operation and the highest level of continuous performance available without the need for helium degassing. Up to four solvent lines may be degassed simultaneously by one unit. One additional channel is reserved for degassing of gradient valve output. Extremely low internal volume of each Systec AFTM channel (480 µl) offers very fast equilibration and very short startup times compared with PTFE degassing channels which have the same degassing efficiency.

Four-way gradient valve is controlled from the pump ECP2010 and 2010H and gradient box ECB2004BP must be ordered with the pump.

Pumps have new sophisticated electronics which allows achieving precise low-pressure gradient.

Up to **6 mobile phase** reservoirs LCT 02 (1L) may be put into the box



7

SPECIFICATION

Degasser	
Degasser channels	5 independent
Degassing process	Gas permeation through a fluoropolymer membrane
Maximum flow rate	10 mL/min.
Maximum operating	0.5 MPa (70 psi)
Degassing capacity	~2 ppm at 1 mL/min.
Degassing efficiency per channel	up to 3 mL/min excellent, up to 5mL/min sufficient
Dead volume of degasser channel	~480 µl per channel for standard channel
Wetted materials	Systec AF [™] , PPS, Glass-filled PTFE, Stainless Steel 316
Gradient valve	
Number of gradient valves	4 (A,B,C,D)
Maximum operating pressure	3 bar (43.5 PSI)
Orifice	0.031" (0.8 mm)
Wetted materials	Simriz [®] , PEEK
Common parameters	
Wetted materials	FEP, PEEK, KEL-F, PTFE, PS, Tefzel™
Connecting thread	1/4"-28 flat bottom
Power supply	100-240V 50/60Hz 40VA
Inside box dimensions (w x h x d)	267 x 56 x 305 mm
Dimensions (w x h x d)	280 x 135 x 498 mm
Main board	3.5" Board Intel Atom Processor E3845 (1.91GHz/4Core
System Memory	SO-DIMM DDR3L-4 GB
Ports	1x USB 2.0; 1 x USB 3.0; 2 x Ethernet 10/100/1000 Mbps; 3x RS232; 1 x RS232/485
Display	1 x VGA; 1 x DDI (DisplayPort, HDMI, DVI)
Storage	SSD 120GB

ECOM spol. s r.o. Trebonicka 239, 252 19 Chrastany u Prahy Czech Republic Tel: +420 221 511 310

www: http://www.ecomsro.com/

Fax: +420 242 498 212

E-mail: info@ecomsro.cz